

pipetman[®]

L



LIMITED EDITION PINK PIPETMAN[®] L

SUPPORTING CANCER RESEARCH WITH EACH PURCHASE



For each Pink PIPETMAN L purchased, \$15 will go to the American Cancer Society

We're Dedicated to Supporting Cancer Research

Worldwide, an estimated 23.6 million people will be diagnosed with cancer by the year 2030.*

Gilson is dedicated to partnering with the scientific community and helping researchers advance their pace of discovering a cure.

With your purchase of the Pink PIPETMAN L, Gilson will donate a portion of the proceeds to the American Cancer Society.

Together we can make a difference in the lives of millions.
Together we can combat cancer.

*Cancer Research UK worldwide cancer statistics



LIMITED EDITION

SHOP YOUR PINK PIPETMAN L AT WWW.GILSON.COM

 **GILSON[®]**



Ultimate Performance


Built for traceability and comfort. PIPETMAN L is designed to increase daily productivity, giving you the ability to pipette longer without fatigue.

Secure Pipetting

Reduce pipetting errors with volume-locking. Avoid accidental volume changes during pipetting cycles, ensuring your accuracy and precision.

Optimized Comfort

Lightweight and light pipetting forces. Right-handed or left-handed? No problem. Adjust the location of the tip ejector to the most comfortable position for you.



TM

Model	Part Number	PIPETMAN DIAMOND Tips		Volume (μL)	Maximum Permissible Errors			
					Gilson		ISO 8655	
					Systematic error (μL)	Random error (μL)	Systematic error (μL)	Random error (μL)
Pink PIPETMAN L P20L	FA10003PINK	D200	DF30	2 10 20	± 0.10 ± 0.10 ± 0.20	≤ 0.030 ≤ 0.050 ≤ 0.060	± 0.2 ± 0.2 ± 0.2	≤ 0.1 ≤ 0.1 ≤ 0.1
Pink PIPETMAN L P200L	FA10005PINK	D200 D300	DF200 DF300	20 100 200	± 0.50 ± 0.80 ± 1.60	≤ 0.20 ≤ 0.25 ≤ 0.30	± 1.6 ± 1.6 ± 1.6	≤ 0.6 ≤ 0.6 ≤ 0.6
Pink PIPETMAN L P1000L	FA10006PINK	D1000 D1200	DF1000 DF1200	100 500 1000	± 3.0 ± 4.0 ± 8.0	≤ 0.6 ≤ 1.0 ≤ 1.5	± 8.0 ± 8.0 ± 8.0	≤ 3.0 ≤ 3.0 ≤ 3.0

Limited edition. Applies to qualifying orders received by Gilson within the limit of available stock.

www.gilson.com/contactus

ISO 9001
CERTIFIED

LT800576/B | ©2019 Gilson
All rights reserved. All product and company names are trademarks of their respective holders. Use of the trademarks in this document does not imply any affiliation with or endorsements by the trademark holders.

GILSON®